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Catalog No.	Product name	volume	원재료
M-CF01-1	Permeabilized solution	15ml	0.2% Tween 20, NaN <sub>3</sub> , pH7.4
M-CF01-2	Blocking solution	15ml	2% 염소 혈청, NaN <sub>3</sub> , pH7.4
M-CF01-3	MARS1 antibody solution	15ml	MARS1 mouse monoclonal Antibody, NaN <sub>3</sub> , BSA, pH7.4
M-CF01-4	Detection antibody solution	15ml	Goat anti-Mouse IgG Secondary Antibody (Alexa Fluor 488), NaN <sub>3</sub> , BSA, pH7.4
M-CF01-5	Mounting Fluid with DAPI	2ml	4',6-diamidino-2-phenylindole, dihydrochloride (DAPI)
M-CF01-6	20X washing solution	50ml * 4ea	20X TBST, 1% Tween 20, pH7.4

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Cell line	# Negative	# Total	Cell line	# Negative	# Total
NIH3T3	3	3	CT26	3	3

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1. Sung Ill Jang et al, New staining method using methionyl-tRNA synthetase 1 antibody for brushing cytology of bile duct cancer. *Gastrointest. Endosc.* 2020; 92:310-9
2. Sunghoon Kim et al, Aminoacyl-tRNA synthetases and tumorigenesis: more than housekeeping. *Nat. Rev. Cancer.* 2011, 11:708-718.
3. Jie Hu et al., Heterogeneity of tumor-induced gene expression changes in the human metabolic network. *Analysis. Nat. Biotechnol.* 2013, 31: 522-529.
4. Eun Young, Kim. et al., Methionyl-tRNA synthetase overexpression is associated with poor clinical outcomes in non-small cell lung cancer. *BMC Cancer.* 2017. 17: 467.
5. Davendra. P.S. Sohail et al., Molecular characteristics of biliary tract cancer. *Crit. Rev. Oncol. Hematol.* 2016, 107: 111-118.
6. Nam Hoon Kwon et al., Stabilization of cyclin-dependent kinase 4 by methionyl-tRNA synthetase in p16INK4a-negative cancer. *ACS Pharmacol. Transl. Sci.* 2018, 1: 21-31.
7. Nam Hoon Kwon et al., Aminoacyl-tRNA synthetases as therapeutic targets. *Nat. Rev. Drug Disc.* 2019, 18: 629-650.
8. Shirley M. Bluethmann et al., Anticipating the "Silver Tsunami": Prevalence trajectories and comorbidity burden among older cancer survivors in the United States. *Cancer Epidemiol. Biomarkers Prev.* 2016, 25:1029-1036.
9. Grand View Research. *Histology and Cytology Market. Market, by Study (Cytology, Histology) and by Cytology (Cervical Cancer, Breast Cancer, Pancreatic Cancer, Bladder Cancer, Bile Duct Cancer), 2017.*
10. Ming-ming Xu, Amrita Sethi. *Diagnosing biliary malignancy. Gastrointest Endosc. Clin. N. Am.* 2015, 25: 677-90.
11. Zaree Babakhanian et al., *Clinics Review Articles. Gastrointestinal Endoscopy Clinics of North America: Endoscopic approach to the patient with biliary tract disease. Editor: Jacques Van Dam, 2013. Vol 23.*
12. In Woong Han et al., Ceruloplasmin as a prognostic marker in patients with bile duct cancer. *Oncotarget.* 2017, 8: 29028-29037.
13. Udayakumar Navaneethan et al., Comparative effectiveness of biliary brush cytology and intraductal biopsy for detection of malignant biliary strictures: a systematic review and meta-analysis. *Gastrointest. Endosc.* 2015, 81: 168-176.

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